



# INSULINK

---

A useful diary for your diabetes

# OUR TEAM

---




**Nicola Dean**

Polimi  
C.P. 10674826



**Marco Fasanella**

Polimi  
C.P. 10617541





## OUR IDEA

---

InsuLink main goal is to give a first support to the patient but only if combined with the doctor supervision.

It is important to underline that this application is only defined by an algorithm, and in this kind of diseases each patient needs ad hoc treatments


# MARKET SIZE

**95%**

Type 2

**5%**

Type 1



Source: WHO

# Functionalities

---

## Scan your food

Scan a given Food BarCode and read all its nutrients!

## Insulin Calculator

An algorithm will retrieve all essential data to calculate the optimal insulin dose

## Calendar


Users can see a well detailed sight of all previous data, just choosing a date from the InsuLink calendar

## Glycemia

Keep track of your daily Glycemia with intuitive charts!

# Charts


## Glycemia Trend




## Macro



## Sport HeatMap




## Sport last week



# Layout

---



## iPad



Dynamic adaption for tablets.  
Currently tested on iPad



## Phone

Compact style with both  
Android and iOS devices  
supported

# Horizontal Layout



## Tablet


iPad and Tablet support also customized horizontal layout



# Demo

## Charts

Useful charts with glycemia Trend and macros



## Shortcuts

Just click a button to reach all the functionalities


## Navigation

Instantly navigate from home menu to the diary or food search

# Search

## Go in detail

Search every kind of food or recipe and get all nutrition details you need



# Scan

## Barcode

Most Famous Food  
is available for scan



## API

Moreover, using  
Premium APIs the

# Glycemia




## Misuration

Save your current glycemia and see its trend

# Insulin Calculator



**Skip**  
Last Glycemia value  
is used



**Last Glycemia**

Insert you last  
misuration

# Architecture



# TEST STRATEGY

---

## Snapshot

Check change in status or pages aspect

## Widget

Test Pages/Components buttons events

## Unit

Test our helper class methods individually

# TEST STRATEGY

---

16



**Test suites**

59



**Normal  
Tests**

33



**Snapshots**



# Testing

---

## Firestore

Check Firestore components are called correctly

## Api

Check api queries

## LocalStorage

Check Local storage function use data correctly

## Redux Reducers

Snapshot reducers behaviour

# Testing

---

## InputChecker

Check Login/Reg errors are thrown correctly

## Pages Snapshots

Snapshots on the page base aspect

## Components

Snapshot custom buttons components

## Pages WidgetTest

Test pages buttons/inputs behaviour

# Future Implementations

---

Some improvements that will make InsuLink even more faster and accurate!



## Machine Learning


Useful to predict  
the correct  
insulin dose

## Premium API

Find even more  
food and recipes

## NFC Technology

Measure your  
glucose directly  
from the App



# THANKS

**Special thanks to Matteo Ciano who designed the app logo. For any business contact you can reach him at: [matteociancio99@gmail.com](mailto:matteociancio99@gmail.com)**

**CREDITS:** This presentation template was created by **Slidesgo**, including icons by **Flaticon**, and infographics & images by **Freepik** and illustrations by **Stories**.

